## Supporting information for

## Molecular level insights into gas oil hydrodenitrogenation by Fourier Transform Ion Cyclotron Resonance Mass Spectrometry

Minh-Tuan Nguyen<sup>1</sup>, Gerhard D. Pirngruber<sup>1\*</sup>, Fabien Chainet <sup>1</sup>, Florian Albrieux<sup>1</sup>, Mélaz Tayakout-Fayolle<sup>2</sup>, Christophe Geantet<sup>3</sup>

<sup>1</sup> IFP Energies nouvelles, Rond-point de l'échangeur de Solaize, BP 3, 69360 Solaize, France

<sup>2</sup> Laboratoire d'Automatique et de Génie des Procédés, Université de Lyon, Université Claude Bernard Lyon 1, CNRS/UCBL, UMR 5007,43 Bd du 11 Novembre 1918, 69622 Villeurbanne Cedex, France

<sup>3</sup> Université de Lyon, IRCELYON, UMR 5256/CNRS-Université Lyon 1, 2 avenue Albert Einstein, F-69626 Villeurbanne, France

## Supporting information 1 : HT-GCxGC coupled by Simulated Distillation analyses of the gas oil feed and the 80% HDN conversion effluent

In order to evaluate the conversion by cracking reactions and hydrocarbon composition after hydrotreatment, the gas oil feed and the hydrotreated effluent at 79.5% HDN conversion were analyzed by high-temperature comprehensive two dimensional chromatography (HT-GCxGC) coupled with simulated distillation. The configuration of GCxGC was described as following:

Parameters	Microfluidic modulator
Instrument	Agilent 7890A
Injector	MMA split mode 400°C split ration 25 $2\mu L$
First column	DB1 HT 15m 0.1mm 0.1µm
Second column	ZB35 HT 5m 0.25mm 0.18µm
Carrier gas	Hydrogène 0.14mL/min et 24mL/min
Oven temperature	100°C 2min 400°C at 2.2°C/min
Modulation time	8.965s
Detector	FID 100Hz

The chromatograms of the 80% HDN conversion effluent is given below:





Retention Time 1 (min)

Figure S2 : FT-ICR/MS raw spectra of the gas oil feed and of the hydrotreated effluents: (a): ESI(+)-FT-ICR/MS spectra, (b): ESI(-)-FT-ICR/MS spectra. From top to bottom: spectra of the gas oil feed, 71%, 74%, 80% and 89% HDN conversion effluents.

(b)



(a)







